



CLINNOVA Gynae 1



The Gynae 1 is the entry-level option to the CLINNOVA Gynae Range. Supplied with a choice of wheel systems, hydraulic or electric height adjustment and classic or premium base, the Gynae 1 can be tailored for a range of specification and budget-driven requirements.

Models

See Page 34-35 for available models



Model NV8463-CLS
Colour Light Grey



Model NV9043
Removable Leg Extension

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted backrest, adjustable from 0° to +90°
- Anti-trap gas assisted backrest to prevent risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46cm to 95cm for individual working heights
- Choice of independent, retractable or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- A pair of multi-adjustable leg supports and arm loops supplied as standard
- Seat section lifts for easy cleaning of framework
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6200** Single foot switch (electric only)
- 6007** Removable head/knee support
- 9043** Removable leg extension
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV8553-PRM
Colour Dark Blue

CLINNOVA Gynae 1 - Models Available



NV8453-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Classic Base



NV8463-CLS

Electric Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Classic Base



NV8553-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Retractable Castors
Classic Base



NV8563-CLS

Electric Height Elevation
Gas Assisted Backrest
Retractable Castors
Classic Base



NV8653-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Central Locking Castors
Classic Base



NV8663-CLS

Electric Height Elevation
Gas Assisted Backrest
Central Locking Castors
Classic Base



NV8453-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Premium Base



NV8463-PRM

Electric Height Elevation
Gas Assisted Backrest
Independently Braking Castors
Premium Base



NV8553-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Retractable Castors
Premium Base



NV8563-PRM

Electric Height Elevation
Gas Assisted Backrest
Retractable Castors
Premium Base



NV8653-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Central Locking Castors
Premium Base



NV8663-PRM

Electric Height Elevation
Gas Assisted Backrest
Central Locking Castors
Premium Base