



CLINNOVA Mobile Hygiene



A selection of easily transportable, variable height changing tables suitable for individuals with additional living needs. Popular within schools for children with special educational needs, but also suitable for adolescents and adults.

Models

See Page 38-41 for available models

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Choice of Small, Medium or Large upholstery surfaces
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Choice of 125mm independent or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6200** Single foot switch (*NV1X61 only*)
- 9002P** Side support rail cushions (Pair)



Model NV1651A-PRM
Colour Lotus Green



Side Support Rails
Fitted as Standard



Model NV1451B-CLS
Colour Dark Blue

CLINNOVA Mobile Hygiene - Models Available



NV1451A-CLS

Hydraulic Height Elevation
Large
Independently Braking Castors
Classic Base



NV1461A-CLS

Electric Height Elevation
Large
Independently Braking Castors
Classic Base



NV1451B-CLS

Hydraulic Height Elevation
Medium
Independently Braking Castors
Classic Base



NV1461B-CLS

Electric Height Elevation
Medium
Independently Braking Castors
Classic Base



NV1451C-CLS

Hydraulic Height Elevation
Small
Independently Braking Castors
Classic Base



NV1461C-CLS

Electric Height Elevation
Small
Independently Braking Castors
Classic Base



NV1451A-PRM

Hydraulic Height Elevation
Large
Independently Braking Castors
Premium Base



NV1461A-PRM

Electric Height Elevation
Large
Independently Braking Castors
Premium Base



NV1451B-PRM

Hydraulic Height Elevation
Medium
Independently Braking Castors
Premium Base



NV1461B-PRM

Electric Height Elevation
Medium
Independently Braking Castors
Premium Base



NV1451C-PRM

Hydraulic Height Elevation
Small
Independently Braking Castors
Premium Base



NV1461C-PRM

Electric Height Elevation
Small
Independently Braking Castors
Premium Base

CLINNOVA Mobile Hygiene - Models Available



NV1651A-CLS

Hydraulic Height Elevation
Large
Central Locking Castors
Classic Base



NV1661A-CLS

Electric Height Elevation
Large
Central Locking Castors
Classic Base



NV1651B-CLS

Hydraulic Height Elevation
Medium
Central Locking Castors
Classic Base



NV1661B-CLS

Electric Height Elevation
Medium
Central Locking Castors
Classic Base



NV1651C-CLS

Hydraulic Height Elevation
Small
Central Locking Castors
Classic Base



NV1661C-CLS

Electric Height Elevation
Small
Central Locking Castors
Classic Base



NV1651A-PRM

Hydraulic Height Elevation
Large
Central Locking Castors
Premium Base



NV1661A-PRM

Electric Height Elevation
Large
Central Locking Castors
Premium Base



NV1651B-PRM

Hydraulic Height Elevation
Medium
Central Locking Castors
Premium Base



NV1661B-PRM

Electric Height Elevation
Medium
Central Locking Castors
Premium Base



NV1651C-PRM

Hydraulic Height Elevation
Small
Central Locking Castors
Premium Base



NV1661C-PRM

Electric Height Elevation
Small
Central Locking Castors
Premium Base