



CLINNOVA Phlebotomy 1



The Phlebotomy 1 is the entry-level option to the CLINNOVA Phlebotomy Range. Supplied with a choice of wheel systems, hydraulic or electric height adjustment and classic or premium base, the Phlebotomy 1 can be tailored for a range of specification and budget-driven requirements.

Models

See Page 58-59 for available models



Model NV9666P-PRM
Colour Canard

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted head and foot sections
- Anti-trap gas assisted sections prevent risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46cm to 95cm for individual working heights
- Choice of Retractable, independent or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Supplied with a pair of multi-adjustable armrest and arm loops
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6200** Single foot switch
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV9566P-CLS
Colour Dark Grey

CLINNOVA Phlebotomy 1 - Models Available



NV9456P-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Independently Braking Castors
Classic Base



NV9466P-CLS

Electric Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Independently Braking Castors
Classic Base



NV9556P-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Retractable Castors
Classic Base



NV9566P-CLS

Electric Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Retractable Castors
Classic Base



NV9656P-CLS

Hydraulic Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Central Locking Castors
Classic Base



NV9666P-CLS

Electric Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Central Locking Castors
Classic Base



NV9456P-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Independently Braking Castors
Premium Base



NV9466P-PRM

Electric Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Independently Braking Castors
Premium Base



NV9556P-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Retractable Castors
Premium Base



NV9566P-PRM

Electric Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Retractable Castors
Premium Base



NV9656P-PRM

Hydraulic Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Central Locking Castors
Premium Base



NV9666P-PRM

Electric Height Elevation
Gas Assisted Backrest
Gas Assisted Foot Section
Central Locking Castors
Premium Base