



## CLINNOVA Phlebotomy 2



The CLINNOVA Phlebotomy 2 is designed for blood sampling procedures, so is fitted with side loops and adjustable upper limb supports to assist with patient positioning and security. The seat section is slightly inclined to encourage the patient to sit back into the chair for additional safety in the event of fainting.

### Models

See Page 62-63 for available models



Model NV9576P-CLS  
Colour Sky Blue

### Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Electric backrest section, gas assisted foot section
- Anti-trap gas assisted sections prevent risk of crushing and harm to patient, user or furniture
- Electric height elevation, raises from low to high in 15 seconds
- Variable height range of 46cm to 95cm for individual working heights
- Choice of Retractable, independent or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Includes a hand switch as standard
- Supplied with a pair of multi-adjustable armrest and arm loops
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

### Optional Accessories

- 6207** Double foot switch
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder  
Fitted as Standard



Model NV9476P-PRM  
Colour Burgundy

## CLINNOVA Phlebotomy 2 - Models Available



### NV9476P-CLS

Electric Height Elevation  
Electric Backrest  
Gas Assisted Foot Section  
Independently Braking Castors  
Classic Base



### NV9576P-CLS

Electric Height Elevation  
Electric Backrest  
Gas Assisted Foot Section  
Retractable Castors  
Classic Base



### NV9676P-CLS

Electric Height Elevation  
Electric Backrest  
Gas Assisted Foot Section  
Central Locking Castors  
Classic Base



### NV9476P-PRM

Electric Height Elevation  
Electric Backrest  
Gas Assisted Foot Section  
Independently Braking Castors  
Premium Base



### NV9576P-PRM

Electric Height Elevation  
Electric Backrest  
Gas Assisted Foot Section  
Retractable Castors  
Premium Base



### NV9676P-PRM

Electric Height Elevation  
Electric Backrest  
Gas Assisted Foot Section  
Central Locking Castors  
Premium Base