



## CLINNOVA Phlebotomy Pro



The CLINNOVA Phlebotomy Pro is intelligently designed for blood sampling procedures, so is fitted with side loops and adjustable upper limb supports to assist with patient positioning and security. The seat section is slightly inclined to encourage the patient to sit back into the chair for additional safety in the event of fainting. Fully electric profiling of height, backrest, footrest and two-way Trendelenburg tilt. The couch also includes a CPR facility and digital display handset, with safe working load of 285kg.

### Models

See Page 66-67 for available models



Model NV9696P-PRM  
Colour Olive

### Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 285Kg (627lbs)
- Fully electric profiling of height, backrest and footrest
- Electric trendelenburg and reverse trendelenburg
- Raises from lowest to highest height in 15 seconds
- Variable height range of 52cm to 100cm to suit individual working heights
- Choice of retractable, independent or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Includes a digital display handset as standard
- Supplied with a pair of multi-adjustable armrest and arm loops
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

### Optional Accessories

- 6004** Breathing hole and plug
- 9210** 10 Button Foot Switch
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder  
Fitted as Standard



Model NV9596P-CLS  
Colour Light Grey

## CLINNOVA Phlebotomy Pro - Models Available



### NV9496P-CLS

Electric Height Elevation & Tilt  
Electric Head Section  
Electric Foot Section  
Independently Braking Castors  
Classic Base



### NV9596P-CLS

Electric Height Elevation & Tilt  
Electric Head Section  
Electric Foot Section  
Retractable Castors  
Classic Base



### NV9696P-CLS

Electric Height Elevation & Tilt  
Electric Head Section  
Electric Foot Section  
Central Locking Castors  
Classic Base



### NV9496P-PRM

Electric Height Elevation & Tilt  
Electric Head Section  
Electric Foot Section  
Independently Braking Castors  
Premium Base



### NV9596P-PRM

Electric Height Elevation & Tilt  
Electric Head Section  
Electric Foot Section  
Retractable Castors  
Premium Base



### NV9696P-PRM

Electric Height Elevation & Tilt  
Electric Head Section  
Electric Foot Section  
Central Locking Castors  
Premium Base