



CLINNOVA Ultra Pro



The CLINNOVA Ultra Pro is intelligently designed for ultrasound, sonography and scanning procedures, offering fully electric profiling of height, backrest, footrest and two-way Trendelenburg tilt. The couch also includes a CPR facility and digital display handset, with a safe working load of 285kg.

Models

See Page 30-31 for available models



Model NV9496-PRM
Colour Burgundy

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 285Kg (627lbs)
- Electric backrest, adjustable from 0° to +90°
- Electric foot section, adjustable from 0° to +90°
- Electric height elevation, raises from 52cm to 100cm in 15 seconds
- Variable height range of 46cm to 95cm for individual working heights
- Choice of independent, retractable or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Includes a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 9210** 10 Button Foot Switch
- 6044** Lithotomy Stirrups
- 6045** Gynaecology Leg Supports
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV9596-CLS
Colour Cream

CLINNOVA Ultra Pro - Models Available



NV9496-CLS

Electric Height Elevation & Tilt
Electric Backrest
Electric Foot Section
Independently Braking Castors
Classic Base



NV9596-CLS

Electric Height Elevation & Tilt
Electric Backrest
Electric Foot Section
Retractable Castors
Classic Base



NV9696-CLS

Electric Height Elevation & Tilt
Electric Backrest
Electric Foot Section
Central Locking Castors
Classic Base



NV9496-PRM

Electric Height Elevation & Tilt
Electric Backrest
Electric Foot Section
Independently Braking Castors
Premium Base



NV9596-PRM

Electric Height Elevation & Tilt
Electric Backrest
Electric Foot Section
Retractable Castors
Premium Base



NV9696-PRM

Electric Height Elevation & Tilt
Electric Backrest
Electric Foot Section
Central Locking Castors
Premium Base