



CLINNOVA Therapy 2 Section (Plus)



A very strong and rigid couch available with a wide range of accessories. The plus head design features two adjustable arm supports for additional positioning techniques. Suitable for supporting very heavy patients.

Models

See Page 14-15 for available models



Model NV2569-PRM
with 9001 Perimeter Foot Switch
Colour Red

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted head section, adjustable from -25° to +45°
- Anti-trap gas assisted head section prevents risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Linked retractable wheels with adjustable foot for uneven floors
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 9001** Perimeter foot switch (NV2569 only)
- 6200** Single foot switch (NV2569 only)
- 9002** Side support rails (Pair)
- 9002P** Side support rail cushions (Pair)
- 9005** Paper Roll Holder
- 6007** Removable head/knee support
- 9008** Integral Face Cradle Cushion



Offset gas strut for unobstructed patient view through breathing hole



Model NV2559-CLS
Colour Olive

CLINNOVA Therapy 2 Section (Plus) - Models Available

NV2559-CLS

Hydraulic Height Elevation
Gas Assisted Head Section
Classic Base



NV2569-CLS

Electric Height Elevation
Gas Assisted Backrest
Hand or Foot Switch
Classic Base



NV2569-CLS-EPFS

Electric Height Elevation
Gas Assisted Head Section
Perimeter Foot Switch
Classic Base



NV2559-PRM

Hydraulic Height Elevation
Gas Assisted Head Section
Premium Base



NV2569-PRM

Electric Height Elevation
Gas Assisted Head Section
Hand or Foot Switch
Premium Base



NV2569-PRM-EPFS

Electric Height Elevation
Electric Head Section
Perimeter Foot Switch
Premium Base

