

Width  
70cm



130cm

65cm

46cm to 95cm



## CLINNOVA Therapy 2 Section (Standard)



A very strong and rigid couch available with a wide range of accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.

### Models

See Page 6-7 for available models



Model NV2567-PRM  
with 9001 Perimeter Foot Switch  
Colour Poterie

### Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted head section, adjustable from -20° to +80°
- Anti-trap gas assisted head section prevents risk of crushing and harm to patient, user or furniture
- Hydraulic pump action or electric height elevation
- Variable height range of 46 to 95cm for individual working heights
- Linked retractable wheels with adjustable foot for uneven floors
- Lock & Go transport mode for easy moving
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Electric models include a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery conceals all visible staples
- Choice of Standard or Standard+ upholstery colours

### Optional Accessories

- 9001** Perimeter foot switch (NV2567 only)
- 6200** Single foot switch (NV2567 only)
- 9002** Side support rails (Pair)
- 9002P** Side support rail cushions (Pair)
- 9005** Paper Roll Holder
- 6007** Removable head/knee support
- 9008** Integral Face Cradle Cushion



Offset gas strut for unobstructed patient view  
through breathing hole



Model NV2557-CLS  
Colour Canard

## CLINNOVA Therapy 2 Section (Standard) - Models Available



### NV2557-CLS

Hydraulic Height Elevation  
Gas Assisted Head Section  
Classic Base



### NV2567-CLS

Electric Height Elevation  
Gas Assisted Backrest  
Hand or Foot Switch  
Classic Base



### NV2567-CLS-EPFS

Electric Height Elevation  
Gas Assisted Head Section  
Perimeter Foot Switch  
Classic Base



### NV2557-PRM

Hydraulic Height Elevation  
Gas Assisted Head Section  
Premium Base



### NV2567-PRM

Electric Height Elevation  
Gas Assisted Head Section  
Hand or Foot Switch  
Premium Base



### NV2567-PRM-EPFS

Electric Height Elevation  
Electric Head Section  
Perimeter Foot Switch  
Premium Base