



## Innovation Deluxe 3 Section



The Innovation Deluxe 3 section couch is designed for use within specialist critical care environments. Electric height, backrest and foot section profiling and 2-way tilt ensure a smooth transition when positioning patients for specialist procedures and treatments. Auto-electric CPR functionality and auto levelling software return the couch quickly to a fully horizontal resuscitation position in the event of an emergency.

### Models

- MG3475** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large independently braked wheel design
- MG3675** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large central locking wheel design with steering facility



Model MG3475 with large independently braked wheel design

### Standard Features

- Electric height range from 47cm to 97cm
- Electric trendelenburg tilt to  $\pm 15^\circ$
- Electric backrest adjustable from  $0^\circ$  to  $+80^\circ$
- Electric footrest adjustable from  $+10^\circ$  to  $-35^\circ$
- Safe working load and lifting capacity of 250Kg (550lbs)
- Auto-electric CPR function via handset
- Auto-intelligent operation
- Vertical lifting design to minimise space and provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps and facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Large TENTE<sup>®</sup>  $\varnothing 125\text{mm}$  castors
- **SEERS ComfortPlus™** Upholstery - features anti-microbial vinyl incorporating pressure relieving design for additional patient comfort with a foam depth of 60mm
- 70cm Standard width upholstery
- Clearance beneath frame for mobile hoist
- Head Cushion & Paper Roll Holder fitted as standard
- Choice of Standard or Standard+ upholstery colours

### Optional Accessories

- 6003** Side support rails (Pair)
- 9002P** Side support rail cushions (Pair)
- 6004** Breathing hole and plug for prone treatments
- 6006** Upper limb support
- 6054** Telescopic IV pole
- 6057** Battery back-up facility



Trendelenburg tilt to  $-15^\circ$



Low patient transfer height



Reverse trendelenburg tilt to  $+15^\circ$

### Available Upholstery Colours:

