



## Innovation Deluxe Dialysis



The Innovation Deluxe Dialysis couch has been designed for use within specialist renal dialysis environments. The couch features auto-electric CPR and fully electric profiling with a 2-way tilt facility to ensure a smooth transition when positioning patients for dialysis procedures and treatments.

### Models

- MG3490** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large independently braked wheel design
- MG3690** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large central locking wheel design with steering facility



Model MG3490 with large independently braked wheel design

### Standard Features

- Electric height range from 47cm to 97cm
- Electric trendelenburg tilt to  $\pm 15^\circ$
- Electric backrest adjustable from  $0^\circ$  to  $+80^\circ$
- Electric footrest adjustable from  $+10^\circ$  to  $-35^\circ$
- Safe working load and lifting capacity of 250Kg (550lbs)
- Auto-electric CPR function via handset
- Auto-intelligent operation
- Vertical lifting design to minimise space and provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps and facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Large TENTE<sup>®</sup>  $\varnothing 125$ mm castors
- **SEERS ComfortPlus™** Upholstery - features anti-microbial vinyl incorporating pressure relieving design for additional patient comfort with a foam depth of 60mm
- 70cm standard width upholstery
- Clearance beneath frame for mobile hoist
- Head Cushion & Paper Roll Holder fitted as standard
- Multi-adjustable arm rests and adjustable upholstered foot board included for patient comfort
- Choice of Standard and Standard+ upholstery colours

### Optional Accessories

- 6054** Telescopic IV pole
- 6057** Battery back-up facility
- 6241** Patient lap tray



Low patient transfer height



2-way trendelenburg tilt facility to  $\pm 15^\circ$



Fully profiled for patient comfort

### Available Upholstery Colours:

