



63	125	70	Standard Dimensions
70*	130	65	*When fitted with 9200/70 Easy Clean Upholstery

## Medicare Mobile 2 Section Couch



Very strong and rigid 2 section mobile treatment couches supplied with side support rails and a paper roll holder. The large wheel system incorporates a steering facility and the couches have excellent height ranges with a low patient transfer height of just 46cm.

### Models

- Independently Braking Castors (IBC)**  
**SM2450** Hydraulic with gas assisted backrest  
**SM2460** Electric with gas assisted backrest  
**SM2470** Electric with electric backrest
- Central Locking Castors (LMWD)**  
**SM2650** Hydraulic with gas assisted backrest  
**SM2660** Electric with gas assisted backrest  
**SM2670** Electric with electric backrest



Model SM2650  
 Shown with optional 6053 Base Cover, 9002P Side Cushions,  
 6054 IV Pole and 6007 Head Cushion  
 Colour Sky Blue

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Gas assisted backrest adjustable from -25° to 75° on SM2450, SM2460, SM2650 and SM2660
- Electric backrest adjustable from -15° to 89° on SM2470/2670
- Gas assisted backrest incorporates anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 46cm to 99cm
- Side Rails & Paper Roll Holder fitted as standard
- Choice of independently or centrally locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 63cm wide Easy Clean Upholstery as standard
- Choice of Standard or Standard+ upholstery colours

### Optional Accessories

- 9002P** Side support cushions (Pair)
- 6004** Breathing hole and plug
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6053** Protective base frame cover
- 6054** Telescopic IV pole
- 6057** Battery back-up facility
- 1WF** Single foot switch for SM2X60
- 2WF** Double foot switch for SM2X70
- 9200/70** 70cm Width Easy Clean Upholstery



Low patient transfer height  
 Shown with optional 6053 Base Cover



Seated position  
 Shown with optional 6053 Base Cover and  
 6006 Upper Limb Support

### Available Upholstery Colours:

