







Medicare Mobile 2 Section Couch







Very strong and rigid 2 section mobile treatment couches supplied with side support rails and a paper roll holder. The large wheel system incorporates a steering facility and the couches have excellent height ranges with a low patient transfer height of just 46cm.

Models

Independently Braking Castors (IBC)

SM2450 Hydraulic with gas assisted backres
SM2460 Electric with gas assisted backrest
SM2470 Electric with electric backrest

Central Locking Castors (LMWD)

SM2650 Hydraulic with gas assisted backrest
SM2660 Electric with gas assisted backrest
SM2670 Electric with electric backrest

Model SM2650

Shown with optional 6053 Base Cover, 9002P Side Cushions, 6054 IV Pole and 6007 Head Cushion

Colour Sky Blue

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Gas assisted backrest adjustable from -25° to 75° on SM2450, SM2460, SM2650 and SM2660
- Electric backrest adjustable from 0° to 85° on SM2470/267
- Gas assisted backrest incorporates anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 46cm to 99cm
- Side Rails & Paper Roll Holder fitted as standard
- Choice of independently or centrally locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 63cm wide Easy Clean Upholstery as standard

Choice of Standard or Standard+ upholstery colours

Optional Accessories

9002P Side support cushions (Pair)

6004 Breathing hole and plug

6006 Upper limb support

6007 Removable head/knee support

6053 Protective base frame cover

6054 Telescopic IV pole

6057 Battery back-up facility

-**1WF** Single foot switch for SM2X60

-2WF Double foot switch for SM2X70

9200/70 70cm Width Easy Clean Upholstery







Seated position Shown with optional 6053 Base Cover an 6006 Upper Limb Support

Available Upholstery Colours:

