





Standard Dimensions





*When fitted with 9300/70 Easy Clean Upholstery

Medicare Mobile 3 Section Couch Low Transfer Height









Very strong and rigid 3 section mobile treatment couches supplied with side support rails and a paper roll holder as standard. The Large Mobile Wheel Design incorporates central locking castors & steering facility, and the couches have excellent height ranges with a low patient transfer height.

Models

Independently Braking Castors (IBC)

SM3450 SM3460 Electric with gas assisted backrest SM3470 Electric with electric backrest

Central Locking Castors (LMWD)

SM3650 Hydraulic with gas assisted backrest SM3660 Electric with gas assisted backrest **SM3670**

> Model SM3650 Colour Sky Blue

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Gas assisted head section adjustable from -25° to +75°
- Choice of hydraulic or electric height elevation

- Side Rails & Paper Roll Holder fitted as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models

Choice of Standard or Standard+ upholstery colours

- Anti-bacterial powder coated frame
- 63cm wide Easy Clean Upholstery as standard
- Telescopic IV pole Choice of independently or centrally locking castors Battery back-up facility
 - Single foot switch for SM3X60 **-2WF** Double foot switch for SM3X70
 - 9300/70 70cm Width Easy Clean Upholstery

Removable head/knee support Protective base frame cover

Optional Accessories

9002P Side support cushions (Pair) **6004** Breathing hole and plug

Upper limb support

Low patient transfer height

Available Upholstery Colours:

