

Medicare Gynaecology Couch




Designed for general gynaecology procedures, these models are available with either hydraulic or electric height operation and feature a wide choice of optional accessories to meet both financial and clinical requirements. All models feature a crescent shaped seat design and are supplied with fully adjustable leg supports, paper roll holder, side support loops and debris tray for gynaecology procedures.

Models

- SM8553** Hydraulic with gas assisted backrest
- SM8563** Electric with gas assisted backrest
- SM8573** Electric height and backrest
- SM8583** Electric height, tilt and backrest



6043
Removable Leg Extension



6051
Fold-Down Leg Extension



Model SM8583
Colour Dark Blue

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric backrest adjustable from 0° to +80°
- Crescent shaped seat design
- Fully adjustable patient leg supports
- Fitted with debris tray, Side Loops & Paper Roll Holder
- Choice of hydraulic or electric height elevation
- Excellent height range from 52cm to 106cm
- Compact base with retractable wheel system
- Adjustable foot for stability on uneven floors
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Hand switch fitted as standard on electric models
- Electric trendelenburg tilt to -15° on model SM8583
- 65cm standard width upholstery
- Choice of Standard and Standard+ upholstery colours

Optional Accessories

- 6007** Head Cushion *(Included on SM8583)*
- 6043** Removable leg extension accessory
- 6051** Fold-down leg extension accessory
- 6052** Protective base frame cover
- 6057** Battery back-up facility
- 1WF** Single foot switch for SM8563
- 2WF** Double foot switch for SM8573
- 3WF** Triple foot switch for SM8583



ST6100
Operators Chair



Model SM8563 Colour light grey

Available Upholstery Colours:

