



## Medicare Multi-Couch



Incorporating both an adjustable angle backrest and a single drop end foot section, these models are suitable for a wide range of seated and supine patient examination procedures. The slightly narrow foot section enables the addition of leg supports or stirrups for gynaecological procedures.

### Models

- SM9556 Hydraulic with gas assisted backrest
- SM9566 Electric with gas assisted backrest
- SM9576 Electric with electric backrest



Model SM9566  
Colour Light Grey



Model 6045 optional leg supports for gynaecology  
Colour may differ from above image



Low patient transfer height



Model 6044 optional lithotomy stirrups

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Gas assisted or electric backrest adjustable from 0° to +80°
- Gas assisted foot section adjustable from 0° to -85°
- Gas assisted sections incorporate an anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 52cm to 106cm
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- 65cm standard width upholstery
- Choice of Standard and Standard+ upholstery colours

### Optional Accessories

- 6004 Breathing hole and plug
- 6005 Paper roll holder
- 6006 Upper limb support
- 6007 Removable head/knee cushion
- 6044 Lithotomy stirrups
- 6045 Leg supports for gynaecology
- 6052 Protective base frame cover
- 6054 Telescopic IV pole
- 6057 Battery back-up facility
- 1WF Single foot switch for SM9566
- 2WF Double foot switch for SM9576

### Available Upholstery Colours:

