


MC6178
45cm to 59cm

MC6179
54cm to 74cm



General Medical Chairs



Suitable for a wide range of medical procedures, these chairs provide comfort and support for clinicians sitting for long periods. The deep padded round seat features a pressure reducing foam membrane and in combination with the fully adjustable backrest, promotes good posture and comfort. There is a choice of standard or high models, with high models also featuring a Foot Support Ring as standard.

Models

MC6178 Round medical chair - standard model
MC6179 Round medical chair - high model
Foot support ring fitted as standard

Features

- Deep padded seat with pressure reducing foam
- Suitable for long periods of sitting
- Weight capacity and safe working load of 150Kg
- Wide choice of models available
- Anti-bacterial and waterproof vinyl upholstery
- Restraint castors fitted as standard
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6178** Anti-static castors
- 6455** Load locking castors
- 6172** Manual locking castors
- 6173** Glides for added stability on high models
- 6175** Height adjustable armrests
- 6174** Foot height control operation
- 6176** Foot support ring for standard model



Model MC6179
Shown with optional foot height control



Model MC6178
Colour Diabolo Cobalt
Colour available at additional cost



Model MC6179
Shown with optional height adjustable armrests
Colour Diabolo Cobalt



Available Upholstery Colours:

