







## Multi Procedures Stool & Chair







Ideally suited for Operating Suites, Ultrasound, Dentistry and Ophthalmic Departments. The 360° swing around arm provides the clinician with an armrest or torso support, using a single movement and offering additional stability during procedures. The ratchet system allows the armrest to be adjusted towards or away from the clinician, providing the user with a variety of options for different procedures. The height of the arm can also be adjusted, ideal if it is being used to support the torso. The backrest models provide additional support and comfort. All models are available with a foot height control as an option.

## Models

MC6160 Chair with backrest - standard model
MC6161 Chair with backrest - high model
Foot support ring fitted as standard

MC6162 Stool - standard model MC6163 Stool - high model

Foot support ring fitted as standard



Fully adjustable 360° swing around arm fitted as standard on all models



## **Features**

- Wide Choice of Models Available
- Intelligent & Ergonomic Design
- 360° Swing around arm fitted on all models and can be used as a backrest or armrest
- Safe Working Load of 150kg
- Upholstery vinyl is Anti-Bacterial, Fireproof to Crib 5 and Waterproof
- Glides (Static Feet) fitted as standard
- Available in a wide range of colours
- Choice of Standard or Standard+ upholstery colours

## **Optional Accessories**

**6178** Anti-static castors

**6455** Load locking castors

**6172** Manual locking castors

173 Glides for added stability on high models

**6174** Foot height control operation

**6176** Foot support ring for standard model





Foot Support Ring fitted as standard on Models MC6161 & MC6163

Model MC6161 Colour Diabolo Cobalt Colour available at additional cost

Available Upholstery Colours:

