


MC6144
45cm to 59cm

MC6145
54cm to 74cm



Ergonomic Saddle Stools



The ergonomic saddle stools have a similar design to the deluxe saddle chairs, designed to encourage positive posture throughout seated procedures. The 2 lever mechanism operation provides independent adjustment for the seat height and tilt. All models feature a polished aluminium five star base and castors.

Models

- MC6144** Saddle stool - standard model
- MC6145** Saddle stool - high model
Foot support ring fitted as standard.




Model MC6145
with Foot Support Ring
Colour Lotus Green



Model MC6144
Colour Lotus Green


Available Upholstery Colours:

STANDARD



Sky Blue — Canard — Dark Blue — Dark Grey — Black

STANDARD+



White — Light Grey — Cream — Olive — Lotus Green — Pepper — Brown — Poterie — Red — Burgundy