









Surgeons & Sonographers Chair









The surgeons and sonographers chairs feature an ergonomically designed seat pad and smaller backrest that will not impede the movement of the arm during scans and procedures. The curvature design and seat pad adjustment encourages pelvic tilt, prevents slouching and promotes positive posture. The chair is standard with a small, compact base to reduce trip hazards and features a wipe clean upholstered underside.

Models

MC6158 MC6159 Surgeons chair - standard model Foot support ring fitted as standard



Fully upholstered backrest and seat pad for ease of cleaning



Features

- Wide choice of models available
- Ergonomically designed seat
- Compact backrest for improved scanning procedures
- Weight capacity and safe working load of 150Kg
- Anti-bacterial and waterproof vinyl upholstery
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- **6178** Anti-static castors
- **6455** Load locking castors
- **6172** Manual locking castors
- Glides for added stability on high models
- Height adjustable armrests
- **6174** Foot height control operation
- **6176** Foot support ring for standard model





Available Upholstery Colours:

