









Sterling 3 Section Couch



Like the 2 Section Model, the Sterling 3 Section Couch is excellent value for money and the 3 section format provides a more versatile and comfortable position for seated procedures and examinations. The backrest support can be raised and fully adjusted from horizontal to +85° and the foot section is adjustable from -20° to +40° and can also be used as the head end for prone treatments with the option of adding a breathing hole and plug. The couch is therefore suitable for a wide range of examination and minor treatment procedures. A hand switch is fitted as standard for height control with a foot switch available as an option.



Standard Features

- Safe working load and lifting capacity of 200Kg
- Gas assisted head/foot section adjustable from -20° to +40°
- Gas assisted backrest adjustable from 0° to +85
- Excellent height range from 46cm to 99cm
- 100mm diameter individual locking castors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- 63cm Wide Easy Clean Upholstery as standard
- Choice of Standard upholstery colours

Optional Accessories

6005 Paper roll holder

6004 Breathing hole and plug

Pair of side support rails

9002P Pair of side support rail cushions

-1WF Single foot switch operation

6007 Matching headrest/ knee support



Available Upholstery Colours:

