

Sterling Changing Table




The Sterling Changing Tables provide a variable height support surface to reduce back strain and posture problems for carers and support staff during routine changing procedures in environments such as special schools, specialist adult centres and paediatric units. These tables are available in two sizes and supplied with side support rails as standard. They are fitted with large wheels with brakes for ease of transport and positioning.

Model

SX1047 Short table electric - 152cm
SX1057 Long table electric - 184cm

Standard Features

- Safe working load and lifting capacity of 200Kg
- Excellent height range from 46cm to 99cm
- Superb low transfer height for patients
- Choice of 2 lengths, 152cm or 184cm
- 100mm diameter individually locking castors
- Side support rails fitted as standard on all models
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on all models
- Standard width upholstery 63cm
- Choice of Standard upholstery colours



Model 9002P Optional Side Rail Cushions can be easily removed from the Side Support Rails for cleaning



Model SX1057
Colour Sky Blue



Model SX1047
Shown with optional 9002P Side Support Rail Cushions
Colour Sky Blue

Available Upholstery Colours:

