







ST7641DL

ST7641

Standard Tilt Tables



These excellent Tilt Tables have been designed to meet the requirements of a wide range of applications for rehabilitation, assessment and examination procedures. Available with the option of "divided leg design" which provides individual harness security for the lower limbs, and options of pommel and leg stabilisers. The standard models include the adjustable work table which is an option on economy models. All models are supplied with individually adjustable foot boards, patient handgrips and harnesses.

Models

ST7641 ST7641DL Standard tilt table

7641DL Standard tilt table (Divided leg)



Model ST7641 Colour Sky Blue



Standard Features

- Safe working load and lifting capacity of 225Kg (495lbs)
- 125mm wheel design, with central locking and steering facility
- Electric height and tilting operation with hand switch control
- Adjustable angle dual foot boards +15° to -30° (4 Positions)
- Height range of 47cm to 101cm
- Tilting angle from 0° to +90°
- Lowers to wheelchair height for ease of patient transfer
- Tilt angle inclinometer fitted as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Supplied with patient handgrips and set of harnesses
- Worktable supplied as standard
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

6057 Battery Back-Up

6058 Hip/Leg Stabilisers (Pair)

6059 Adjustable Pommel Support*

6007 Head Support Cushion

6062 Standard Harness Set6063 Dual Leg Harness Set*

6064 Pair Dual Knee Harnesses*

6065 Standard Harness (Wide)

6066 Standard Harness (Narrow)

* (will fit Divided Leg version only)



Model ST7641 Colour Sky Blue



Available Upholstery Colours:

