





Therapy 2 Section Plus Head Couch



The plus head design is ideal for neck mobilisation as well as manual traction techniques. These couches feature a short head section with height adjustable armrests for prone treatments, promoting relaxed shoulders and a more comfortable patient. The couch has an excellent height range of 45cm to 98cm and is fitted with a breathing hole and plug as standard.

Models

ST2559 Hydra ST2569 Electr

Hydraulic with gas assisted head section Electric with gas assisted head section



Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- The gas assisted head section incorporates anti-trap design gas strut
- Gas assisted head section adjustable from -25° to +40
- Adjustable height armrests as standard
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hois
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

-EPFS Perimeter Foot Switch

-1WF Single Foot Switch

6006 Upper Limb Support

5002 Side Support Rails

9002P Side Support Rail Cushions

6005 Paper Roll Holder

6007 Removable Head/Knee Support





Low patient transfer height

Available Upholstery Colours:

