



Therapy 2 Section Standard Head Couch



A very strong and rigid 2 section couch with a wide range of optional accessories available and an extensive choice of colours. This couch is suitable for supporting very heavy patients for examinations, minor procedures and physiotherapy treatments. This 2 section couch is available as a hydraulic or electric model, both models feature a gas assisted backrest that can be adjusted from -25° to +75° and has an excellent height range from 45cm to 98cm.

Models

- ST2557 Hydraulic with gas assisted head section
- ST2567 Electric with gas assisted head section



Model ST2567
Shown with optional Foot Switch
Colour Light Grey

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- The gas assisted head section incorporates anti-trap design gas strut
- Gas assisted head section adjustable from -25° to +75°
- Choice of hydraulic or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- 63cm wide Easy Clean Upholstery as standard
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- EPFS Perimeter Foot Switch
- 1WF Single Foot Switch
- 6006 Upper Limb Support
- 6002 Side Support Rails
- 9002P Side Support Rail Cushions
- 9005 Paper Roll Holder
- 6007 Removable Head/Knee Support



Negative head section to -25°



Low patient transfer height

Available Upholstery Colours:

