





ST4552

ST4562

ST4542

125 135

> ST4552W ST4562W ST4542W

Height range 46cm to 100cm

Therapy Bobath Couch



By tradition, Bobath treatment techniques require a suitably wide treatment surface with a safe and capable lifting capacity for both patient and therapist. The Bobath Therapy range of couches has more than achieved this with a choice of 2 sizes and higher than standard safe working load of 250Kg. See deluxe model for 325Kg safe working load.

Models ST4552 ST4562

ST4542



Standard Features

- Safe working load and lifting capacity of 250Kg (550lbs)
- The gas assisted backrest section incorporates anti-trap design gas strut
- Choice of hydraulic or electric height elevation
- on uneven floors and surfaces
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

Breathing Hole & Plug Battery Back-Up

Single Foot Switch for ST4562/ST4562W Double Foot Switch for ST4542/ST4542W

Perimeter Foot Switch Side Support Rails

9002P Side Support Rail Cushions

Low patient transfer height

Available Upholstery Colours:

