



Therapy Mat Table



The Mat Therapy Tables are very similar to the Bobath models in terms of specification, but without the backrest facility. The single section top design provides a large and stable variable height support surface for a wide range of exercises and Bobath rehabilitation procedures.

Models

ST5551 Hydraulic, 105cm Wide ST5561 Electric, 105cm Wide ST5551W Hydraulic, 125cm Wide ST5561W Electric, 125cm Wide ST5561L Electric, 150cm Wide





Standard Features

- Safe working load and lifting capacity of 250Kg (550lbs)
- Choice of widths available 105cm, 125cm or 150cm
- Choice of hydraulic or electric height elevation on 105cm and 125cm models
- Excellent electric neight range from 46cm to 100cm
- Retractable wheel design and adjustable foot for stability on uneven floors
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard
- Bariatric models with a safe working load of 325Kg are available upon request
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

6057 Battery Back-Up
-1WF Single Foot Switch
-EPFS Perimeter Foot Switch
6074 Side Support Rails
9002P Side Support Rail Cushions



Available
Upholstery
Colours:

26

