





Therapy Traction Table



The traction tables have been specially designed for lumbar and cervical traction. They feature an integral rolling top design to facilitate effective traction procedures. The tables are fitted with a machine mount as standard.

Models

ST6557 Hydraulic traction table with machine mount ST6567



Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Integral rolling top design to provide effective application of traction forces

- Head section includes breathing hole and plug
- Gas assisted backrest and foot sections incorporate anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 50cm to 101cm
- Fixation point for optional cervical traction accessory

- Standard width 70cm upholstery
- Choice of Standard or Standard+ upholstery colours



-1WF Single Foot Switch

Cervical In Sitting Accessory

Pelvic Lumbar Traction Harness

Thoracic Harness

Thoracic Anchor Straps

Cervical Traction Head Halter

Cervical Pillow

Flexion Stool

Spreader Bar





Available Upholstery Colours:

