


## Tilt Table Pro with Emergency Override



These intelligent tilt tables have been designed & developed with medical professionals to meet requirements of syncope testing. The Emergency Override allows the clinician to return the patient to horizontal immediately, for situations such as fainting or cardiac arrest. A digital display handset is included as standard, giving live feedback of the height and angle that the table is adjusted to. The tables are available in either single or divided leg design.

### Models

- ST7645 Tilt Table Pro - Single Leg
- ST7645DL Tilt Table Pro - Divided Leg



Emergency Override fitted on all models



Model ST7645  
Colour Sky Blue



Lowers to 49cm for ideal transfer height



Set of harnesses & worktable included as standard

### Standard Features

- Safe working load and lifting capacity of 225kg (495lbs)
- Electric height elevation from 49cm to 104cm
- Electric tilt from 0° to 90° in only 15 seconds
- Emergency Down Facility to return patient to horizontal
- Harness set, handgrips and worktable included as standard
- 70cm wide Easy Clean Upholstery fitted as standard
- Digital Display Handset shows live feedback of height (cm) and tilt angle of the table
- Memory function with 3 programmable positions
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Large Mobile Wheel Design with central locking castors
- Anti-bacterial powder coated frame
- Full length accessory fitting bar for optional lateral supports
- Lifetime frame warranty & 5 year warranty on moving parts
- Choice of Standard or Standard+ upholstery colours

### Optional Accessories

- 6057 Battery Back-Up
- 6058 Lateral Supports for Stabilisation
- 6059 Adjustable Pommel Support\*
- 6007 Head Support Cushion
- 6062 Single Leg Harness Set
- 6063 Divided Leg Harness Set\*
- 6064 Pair of Dual Knee Harnesses\*
- 6065 Standard Harness (Wide)
- 6066 Standard Harness (Narrow)

\* for ST7645DL only

### Available Upholstery Colours:

